

PENTAPHARM

Pefabloc[®]TH (α -NAPAP)

Description: The benzamidine moiety is a key structure for the development of inhibitors of trypsin-like enzymes. It mimics the protonated side-chain of the basic amino acids Arg and Lys which is hydrolysed by these enzymes. Pefabloc[®]TH (NAPAP) is one of the most potent and selective competitive inhibitors of thrombin.

Application: Pefabloc[®]TH can be used in diagnostic systems, analytical applications, research and industrial purification processes to exclude undesired thrombin activity. It can also be used as a potent anticoagulant in *in vitro* test systems.

Formula: C₂₇H₃₁O₄N₅S·AcOH **MW:** 581.7

Chem. name: N α -(2-Naphthylsulfonylglycyl)-4-amidino-(D,L)-phenylalanine piperidide acetate (NAPAP)

Solubility: 1.5 mg/ml (2.6 mM) in H₂O

Storage: May be used by the expiry date given on the label when stored unopened, protected from moisture, in the dark, 2 - 8°C. Avoid contamination of the reagents by micro-organisms. Shipment of product does not require cooling during the time of transportation.

K_i constants: K_i values for the inhibition of several trypsin-like serine proteases by Pefabloc[®]TH

Enzyme	K _i [μ M]
Trypsin	0.69
Thrombin (bovine)	0.006
Factor Xa	7.9
Factor XIIa	450
Plasmin	30
uPA	230
sc-tPA	430
tc-tPA	70
Plasma kallikrein	14.4
Glandular kallikrein	93
Acrosin	2.9
Batroxobin (B. atrox)	1.7
Batroxobin (B. moojeni)	3.8
PCa	4.8
Tryptase	45
Chymotrypsin	625

References: Stürzebecher J, Walsmann P, Voigt B, Wagner G. Inhibition of bovine and human thrombins by derivatives of benzamidine. *Thromb Res* 1984; 36: 457-65.

Package size: Bulk [mg] / minimum order 5mg

Code: 381-01

FOR RESEARCH USE ONLY. NOT FOR HUMAN USE OR DRUG USE.